

# Inhalation Therapy Cart

## Model ITC12C


SUPER ERECTA SHELF®

### Inhalation Therapy Cart

- Super Erecta Shelf® Construction:**  
Allows custom height adjustment and positive locking without the use of nuts, bolts or setscrews and without the use of special tools.  
**Note:** To maintain the conductivity of the cart, aluminum split sleeves must be used rather than plastic.
- PVC-Coated Top & Bottom Inserts:** The compartmented inserts, which snap over the top and bottom frames, have a resilient PVC coating which aids in reducing vibration and noise.
- Chrome-Plated Finish** on handles, frames and shelf.
- Accommodates Any Combination** of twelve D & E cylinders. Each PVC-coated insert has twelve compartments, each measuring 4½" (114mm) square.
- 5" (127mm) Conductive Casters.**  
See caster NOTICE on assembly instruction sheet.



ITC12C



Copyright © 1984 InterMetro Industries Corp.

TIM CAMPBELL 800-400-7500

# ITC12C Cart

## DIMENSIONS

| Cat. No. | Width                          |      | Length |      | Height |      | Approx. Pkd. Wt.               |      |
|----------|--------------------------------|------|--------|------|--------|------|--------------------------------|------|
|          | (in.)                          | (mm) | (in.)  | (mm) | (in.)  | (mm) | (lbs.)                         | (kg) |
| ITC12C   | 20 <sup>1</sup> / <sub>2</sub> | 520  | 26     | 660  | 39     | 990  | 40 <sup>1</sup> / <sub>2</sub> | 18   |

TIM CAMPBELL 800-400-7500