

Inhalation Therapy Cart

Model ITC12C


SUPER ERECTA SHELF®

Inhalation Therapy Cart

- Super Erecta Shelf® Construction:**
Allows custom height adjustment and positive locking without the use of nuts, bolts or setscrews and without the use of special tools.
Note: To maintain the conductivity of the cart, aluminum split sleeves must be used rather than plastic.
- PVC-Coated Top & Bottom Inserts:** The compartmented inserts, which snap over the top and bottom frames, have a resilient PVC coating which aids in reducing vibration and noise.
- Chrome-Plated Finish** on handles, frames and shelf.
- Accommodates Any Combination** of twelve D & E cylinders. Each PVC-coated insert has twelve compartments, each measuring 4½" (114mm) square.
- 5" (127mm) Conductive Casters.**
See caster NOTICE on assembly instruction sheet.



ITC12C



Copyright © 1984 InterMetro Industries Corp.

TIM CAMPBELL 800-400-7500

ITC12C Cart

DIMENSIONS

Cat. No.	Width		Length		Height		Approx. Pkd. Wt.	
	(in.)	(mm)	(in.)	(mm)	(in.)	(mm)	(lbs.)	(kg)
ITC12C	20 ¹ / ₂	520	26	660	39	990	40 ¹ / ₂	18

TIM CAMPBELL 800-400-7500